



2000 Multipurpose Patient Chair

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

2000 Exam Chair for Specialty Procedures

Ergonomic Features

Handrails and footrest puts patients at ease and increase safety.

Square base provides crucial stability during procedures.

Upholstered back adjusts vertically for patient comfort.

Power Base with Corded Footswitch

Corded footswitch raises and lowers the 2000 chair, and can be placed anywhere around the base.

Rotation Lock

The versatile 2000 chair swivels 360° and locks in place with a rotation lever.

Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Black upholstery is standard, or choose from 14 standard vinyl colors and endless custom selections - call Reliance at 1-800-735-0357 for a color card.

UL Mark of Quality

The 2000 has passed all rigorous Laboratories tests, achieving CSA, C/US certificate of compliance to assure safety and durability.

Tech Specs

Base: 18" square. Seat height: 19 ¾" lowered, 28 ¾" raised. Shipping weight: 98 lbs. Electrical: 115 volts, 60 Hz, 2.4 amperes.

